


## p38 MAPK Antibody / MAPK14 [clone 5A1-C5-F11] (F54021)

| Catalog No.  | Formulation                                    | Size   |
|--------------|------------------------------------------------|--------|
| F54021-0.1ML | In PBS with 50% glycerol and 0.03% ProClin 300 | 0.1 ml |

**Bulk quote request**

|                           |                                                            |
|---------------------------|------------------------------------------------------------|
| <b>Availability</b>       | 1-3 business days                                          |
| <b>Species Reactivity</b> | Human, Mouse, Rat, Primate                                 |
| <b>Format</b>             | Purified                                                   |
| <b>Clonality</b>          | Monoclonal (mouse origin)                                  |
| <b>Isotype</b>            | Mouse IgG1                                                 |
| <b>Clone Name</b>         | 5A1-C5-F11                                                 |
| <b>Purity</b>             | Protein G affinity                                         |
| <b>UniProt</b>            | Q16539                                                     |
| <b>Applications</b>       | Western Blot : 1:500                                       |
| <b>Limitations</b>        | This p38 MAPK antibody is available for research use only. |



Western blot testing of mouse C2C12, rat PC12, human HepG2, monkey COS7, rat C6, mouse NIH3T3 and human HeLa cell lysates with p38 MAPK antibody at 1:500. Predicted molecular weight: 38-41 kDa.

### Description

Serine/threonine kinase which acts as an essential component of the MAP kinase signal transduction pathway. MAPK14 is one of the four p38 MAPKs which play an important role in the cascades of cellular responses evoked by extracellular stimuli such as proinflammatory cytokines or physical stress leading to direct activation of transcription factors. Accordingly, p38 MAPKs phosphorylate a broad range of proteins and it has been estimated that they may have approximately 200 to 300 substrates each.

### Application Notes

The stated application concentrations are suggested starting points. Titration of the p38 MAPK antibody may be required

due to differences in protocols and secondary/substrate sensitivity.

## **Immunogen**

A human recombinant partial protein was used as the immunogen for this MAPK14 / p38 MAPK antibody.

## **Storage**

Store the p38 MAPK antibody at -20oC.