


p63 Antibody (R35693)


Catalog No.	Formulation	Size
R35693-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8626
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This p63 antibody is available for research use only.



Western blot testing of human heart lysate with p63 antibody at 0.3ug/ml. Expected molecular weight: ~77kDa.



Western blot of mouse (A) and rat (B) heart lysate with p63 antibody at 0.3ug/ml. Expected molecular weight: ~77kDa.

Description

Additional name(s) for this target protein: Tumor protein p63; TP63

Application Notes

Optimal dilution of the p63 antibody should be determined by the researcher.

1. This p63 antibody will detect isoform 1/4.

Immunogen

Amino acids RNKQQRKEEGE were used as the immunogen for this p63 antibody.

Storage

Aliquot and store the p63 antibody at -20oC.