


p75NTR Antibody / NGF Receptor / CD271 [clone DBO-14] (RQ5203)


Catalog No.	Formulation	Size
RQ5203	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

Bulk quote request

Availability	1-2 weeks
Species Reactivity	Human, Rat
Format	Purified
Clonality	Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	DBO-14
Purity	Affinity purified
UniProt	P08138
Applications	Western Blot : 1:10000-1:20000 Immunohistochemistry (FFPE) : 1:50-1:200
Limitations	This p75NTR antibody is available for research use only.



Western blot testing of rat 1) C6 and 2) PC-12 cell lysate with p75NTR antibody. Expected molecular weight: 45-75 kDa depending on glycosylation level.



IHC staining of FFPE human glioma with p75NTR antibody. HIER: boil tissue sections in pH6, 10mM citrate buffer, for 10-20 min and allow to cool before testing.

Application Notes

Optimal dilution of the p75NTR antibody should be determined by the researcher.

Immunogen

A synthetic peptide specific to human NGF Receptor / p75NTR / CD271 was used as the immunogen for the p75NTR antibody.

Storage

Store the p75NTR antibody at -20oC.