


PAK1 Antibody (R34929)


Catalog No.	Formulation	Size
R34929-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5058
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This PAK1 antibody is available for research use only.



Western blot testing of rat brain lysate with PAK1 antibody at 1ug/ml. Expected molecular weight 60~70 kDa.



Western blot testing of human brain (cerebral cortex) lysate with PAK1 antibody at 0.5 ug/ml. Expected molecular weight 60~70 kDa.

Description

Additional name(s) for this target protein: P21/Cdc42/Rac1-activated kinase 1

Application Notes

Optimal dilution of the PAK1 antibody should be determined by the researcher.

Immunogen

Amino acids NTEKQKKKPKMSDE were used as the immunogen for this PAK1 antibody.

Storage

Aliquot and store the PAK1 antibody at -20oC.