


## PAPD5 Antibody / TRF4-2 (R36141)

Catalog No.	Formulation	Size
R36141-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
UniProt	Q68ED3
Gene ID	64282
Applications	Western Blot : 0.2-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This PAPD5 antibody is available for research use only.



Western blot testing of mouse spleen lysate with PAPD5 antibody at 0.2ug/ml.  
Predicted/observed molecular weight (mouse): ~70kDa.

## Description

Additional name(s) for this target protein: PAP associated domain containing 5; TRF4-2

## Application Notes

Optimal dilution of the PAPD5 antibody should be determined by the researcher.

## Immunogen

Amino acids TDEVATYRDWISKQ were used as the immunogen for this PAPD5 antibody.

## Storage

Aliquot and store the PAPD5 antibody at -20oC.