


PARP3 Antibody (F40077)

Catalog No.	Formulation	Size
F40077-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40077-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	Q8CFB8
Applications	Western Blot : 1:1000
Limitations	This PARP3 antibody is available for research use only.



Description

PARP3, poly(ADP-ribosyl)transferase 3, belongs to the PARP family. These enzymes modify nuclear proteins by poly-ADP-ribosylation, which is required for DNA repair, regulation of apoptosis, and maintenance of genomic stability. This protein is preferentially localized to the daughter centriole throughout the cell cycle.

Application Notes

Titration of the PARP3 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 451-481 from the mouse protein was used as the immunogen for this PARP3 antibody.

Storage

Aliquot the PARP3 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.