


PAX8 Antibody / Paired box 8 [clone RM436] (R20451)

Catalog No.	Formulation	Size
R20451-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM436
Purity	Protein A purified from animal origin-free supernatant
UniProt	Q06710
Localization	Nuclear
Applications	Immunohistochemistry (FFPE) : 1:100 -1:400
Limitations	This recombinant PAX8 antibody is available for research use only.



IHC staining of FFPE human thyroid tissue with recombinant PAX8 antibody at 1:400.



IHC staining of FFPE human kidney tissue with recombinant PAX8 antibody at 1:400.

Description

This antibody reacts to PAX8.

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant PAX8 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A peptide corresponding to the C-terminus of PAX8 was used as the immunogen for the recombinant PAX8 antibody.

Storage

Store the recombinant PAX8 antibody at -20oC.