


## PCDH17 Antibody (R35803)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R35803-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                      |                                                             |
|----------------------|-------------------------------------------------------------|
| Availability         | 1-3 business days                                           |
| Species Reactivity   | Human, Mouse                                                |
| Predicted Reactivity | Cow, Pig, Rat                                               |
| Format               | Antigen affinity purified                                   |
| Clonality            | Polyclonal (goat origin)                                    |
| Isotype              | Goat Ig                                                     |
| Purity               | Antigen affinity                                            |
| Gene ID              | 27253                                                       |
| Applications         | Western Blot : 0.1-0.3ug/ml<br>ELISA (peptide) LOD : 1:8000 |
| Limitations          | This PCDH17 antibody is available for research use only.    |



Western blot testing of fetal mouse brain lysate with PCDH17 antibody at 0.2ug/ml.  
Predicted molecular weight: ~126kDa, observed here at ~150kDa.

## Description

Additional name(s) for this target protein: Protocadherin 17

## Application Notes

Optimal dilution of the PCDH17 antibody should be determined by the researcher.

## Immunogen

Amino acids RVFADVHSRASRDS were used as the immunogen for this PCDH17 antibody.

## Storage

Aliquot and store the PCDH17 antibody at -20oC.