


PCSK9 Antibody (R36265)

Catalog No.	Formulation	Size
R36265-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	255738
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This PCSK9 antibody is available for research use only.



Western blot of HEK293 lysate overexpressing human PCSK9 with MYC tag probed with PCSK9 antibody [1ug/ml] in Lane A and probed with anti-MYC tag [1/1000] in lane C. Mock-transfected HEK293 probed with PCSK9 antibody [1ug/ml] in Lane B.

Description

Additional name(s) for this target protein: Proprotein convertase subtilisin/kexin type 9

Application Notes

Optimal dilution of the PCSK9 antibody should be determined by the researcher.

Immunogen

Amino acids TRFHRQASKCDSHGT were used as the immunogen for this PCSK9 antibody.

Storage

Aliquot and store the PCSK9 antibody at -20oC.