


PD1 Antibody (V3384)

Catalog No.	Formulation	Size
V3384-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	100 ug
V3384-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	20 ug
V3384SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Protein A affinity chromatography
UniProt	Q15116
Localization	Cell surface & cytoplasmic
Applications	Immunohistochemistry (FFPE) : 1-2ug/ml for 30 min at RT
Limitations	This PD1 antibody is available for research use only.



IHC testing of FFPE human tonsil with PD1 antibody. Required HIER: steam sections in 10mM Tris with 1mM EDTA, pH 9, for 10-20 min.

Description

PDCD1 (Programmed cell death-1 protein), also called PD1 and CD279, is a type I transmembrane receptor and a member of the immunoglobulin gene superfamily. It is expressed on activated T-cells, B-cells, and myeloid cells. Anti-PDCD1 is a marker of angioimmunoblastic lymphoma and suggests a unique cell of origin for this neoplasm. Unlike CD10

and BCL6, PDCD1 is expressed by few B-cells, so may be a more specific and useful diagnostic marker in angioimmunoblastic lymphoma.

Application Notes

Titering of the PD1 antibody may be required for optimal performance.

Immunogen

A full-length human protein was used as the immunogen for the PD1 antibody.

Storage

Store the PD1 antibody at 2-8oC (with azide) or aliquot and store at -20oC or colder (without azide).