


PEDF Antibody / Pigment epithelium-derived factor / SERPINF1 (R34488)

Catalog No.	Formulation	Size
R34488-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5176
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This PEDF antibody is available for research use only.



Western blot testing of human albumin-depleted plasma [A] and serum [B] lysates with PEDF antibody at 0.3ug/ml. Predicted molecular weight: ~46kDa.

Description

Additional name(s) for this target protein: Pigment epithelium-derived factor; SERPINF1

Application Notes

Optimal dilution of the PEDF antibody should be determined by the researcher.

Immunogen

Amino acids KITGKPIKLTQVEHR were used as the immunogen for this PEDF antibody.

Storage

Aliquot and store the PEDF antibody at -20oC.