


Peripherin Antibody (R35323)

Catalog No.	Formulation	Size
R35323-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	19132
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Peripherin antibody is available for research use only.



Western blot testing of rat colon lysate with Peripherin antibody at 0.1ug/ml. Expected molecular weight ~55kDa. Both bands are blocked by the immunizing peptide.

Description

Additional name(s) for this target protein: Peripherin1; Prph

Application Notes

Optimal dilution of the Peripherin antibody should be determined by the researcher.

Immunogen

Amino acids QKEQHSDLKSSIH were used as the immunogen for this Peripherin antibody.

Storage

Aliquot and store the Peripherin antibody at -20oC.