


Peroxiredoxin 2 Antibody (R34939)

Catalog No.	Formulation	Size
R34939-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7001
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This Peroxiredoxin 2 antibody is available for research use only.



Western blot testing of human brain lysate with Peroxiredoxin 2 antibody at 1ug/ml.
Predicted molecular weight: ~22kDa.

Description

Additional name(s) for this target protein: PRDX2

Application Notes

Optimal dilution of the Peroxiredoxin 2 antibody should be determined by the researcher.

1. This Peroxiredoxin 2 antibody will detect isoforms a and b.

Immunogen

Amino acids NVDDSKEYFSKHN were used as the immunogen for this Peroxiredoxin 2 antibody.

Storage

Aliquot and store the Peroxiredoxin 2 antibody at -20oC.