


## Peroxiredoxin 6 Antibody (R33709)


| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R33709-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                    |                                                                   |
|--------------------|-------------------------------------------------------------------|
| Availability       | 1-3 business days                                                 |
| Species Reactivity | Human, Mouse                                                      |
| Format             | Antigen affinity purified                                         |
| Clonality          | Polyclonal (goat origin)                                          |
| Isotype            | Goat Ig                                                           |
| Purity             | Antigen affinity                                                  |
| Gene ID            | 9588                                                              |
| Applications       | Western Blot : 0.1-0.3ug/ml<br>ELISA (peptide) LOD : 1:128000     |
| Limitations        | This Peroxiredoxin 6 antibody is available for research use only. |



Western blot testing of human liver [A] and heart [B] lysates with Peroxiredoxin 6 antibody at 0.1ug/ml. Expected molecular weight: ~25kDa.



Western blot testing of mouse liver lysate with Peroxiredoxin 6 antibody at 0.1ug/ml. Expected molecular weight: ~25kDa.

## Description

Additional name(s) for this target protein: PRDX6

## **Application Notes**

Optimal dilution of the Peroxiredoxin 6 antibody should be determined by the researcher.

## **Immunogen**

Amino acids EANTTVGRIRFHD were used as the immunogen for this Peroxiredoxin 6 antibody.

## **Storage**

Aliquot and store the Peroxiredoxin 6 antibody at -20oC.