


## PGA4 Antibody (F40025)


Catalog No.	Formulation	Size
F40025-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40025-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

**Bulk quote request**

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Predicted Reactivity</b>	Rabbit
<b>Format</b>	Antigen affinity purified
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	P0DJD7
<b>Applications</b>	Western Blot : 1:1000 IHC (Paraffin) : 1:10-1:50
<b>Limitations</b>	This PGA4 antibody is available for research use only.



PGA4 antibody western blot analysis in WiDr lysate.



PGA4 antibody immunohistochemistry analysis in formalin fixed and paraffin embedded human stomach tissue.

## Description

PGA4 shows particularly broad specificity; although bonds involving phenylalanine and leucine are preferred, many others are also cleaved to some extent.

## Application Notes

Titration of the PGA4 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

## Immunogen

A portion of amino acids 353-381 from the human protein was used as the immunogen for this PGA4 antibody.

## Storage

Aliquot the PGA4 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.