


## PGP9.5 Antibody (R35289)


Catalog No.	Formulation	Size
R35289-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	7345
Applications	Western Blot : 0.01-0.03ug/ml Enzyme ImmunoAssay : binds recombinant protein ELISA (peptide) LOD : 1:16000
Limitations	This PGP9.5 antibody is available for research use only.



Western blot testing of human hippocampus lysate with PGP9.5 antibody at 0.01ug/ml.  
Predicted molecular weight ~25kDa.



Western blot of mouse (A) and rat (B) brain lysates with PGP9.5 antibody at 0.01ug/ml.  
Predicted molecular weight ~25kDa.

## Description

Additional name(s) for this target protein: PARK5, Ubiquitin carboxyl-terminal esterase L1, UCHL1

## Application Notes

Optimal dilution of the PGP9.5 antibody should be determined by the researcher.

## Immunogen

Amino acids QHENFRKKQIE were used as the immunogen for this PGP9.5 antibody.

## Storage

Aliquot and store the PGP9.5 antibody at -20oC.