


Phospho-Src Antibody (pY530) (F48453)


Catalog No.	Formulation	Size
F48453-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F48453-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Chicken, Mouse, Rat, Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P12931
Applications	Dot Blot : 1:500 IHC (Paraffin) : 1:50-1:100
Limitations	This phospho-Src antibody is available for research use only.



Dot blot analysis of phospho-Src antibody. 50ng of phos-peptide or nonphos-peptide per dot were spotted.



IHC analysis of FFPE human hepatocarcinoma tissue stained with the phospho-Src antibody.

Description

This protein has tyrosine kinase activity and belongs to the src family of proteins.

Application Notes

Titration of the phospho-Src antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

This phospho-Src antibody was produced from rabbits immunized with a KLH conjugated synthetic phosphopeptide corresponding to amino acid residues surrounding pY530 of human Src.

Storage

Aliquot the phospho-Src antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.