

PIP5K1B Antibody (F40169)


Catalog No.	Formulation	Size
F40169-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F40169-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Purified
UniProt	O14986
Applications	Western Blot : 1:1000 IHC (Paraffin) : 1:50-1:100
Limitations	This PIP5K1B antibody is available for research use only.

250
150
100
75
50
37
25
20
15

Western blot analysis of PIP5K1B antibody and human placenta cell lysate



IHC analysis of FFPE human brain tissue stained with PIP5K1B antibody

Description

Protein kinases are enzymes that transfer a phosphate group from a phosphate donor, generally the γ phosphate of ATP, onto an acceptor amino acid in a substrate protein. By this basic mechanism, protein kinases mediate most of the signal transduction in eukaryotic cells, regulating cellular metabolism, transcription, cell cycle progression, cytoskeletal rearrangement and cell movement, apoptosis, and differentiation. With more than 500 gene products, the protein kinase family is one of the largest families of proteins in eukaryotes. The family has been classified in 8 major groups based on sequence comparison of their tyrosine (PTK) or serine/threonine (STK) kinase catalytic domains.

Application Notes

Titration of the PIP5K1B antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 1-30 from the human protein was used as the immunogen for this PIP5K1B antibody.

Storage

Aliquot the PIP5K1B antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.