

Pirin Antibody (R34139)


Catalog No.	Formulation	Size
R34139-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	8544
Applications	Western Blot : 0.3-1.0ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This Pirin antibody is available for research use only.



Western blot testing of human heart lysate with Pirin antibody at 0.3ug/ml.
Predicted/observed molecular weight: ~32kDa.



Western blot of HEK293 lysate overexpressing PIR probed with Pirin antibody (mock transfection in lane 1). Predicted/observed molecular weight: ~32kDa.

Description

Additional name(s) for this target protein: Iron-binding nuclear protein; PIR

Application Notes

Optimal dilution of the Pirin antibody should be determined by the researcher.

Immunogen

Amino acids GSSKKVTLSVLSREC were used as the immunogen for this Pirin antibody.

Storage

Aliquot and store the Pirin antibody at -20oC.