


PITPNB Antibody (R35731)

Catalog No.	Formulation	Size
R35731-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Cow, Dog, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5306
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This PITPNB antibody is available for research use only.



Western blot testing of 293 lysate with PITPNB antibody at 1ug/ml. Predicted/observed molecular weight: ~32kDa.

Description

Additional name(s) for this target protein: PITPN, PITP alpha, phosphatidylinositol transfer protein

Application Notes

Optimal dilution of the PITPNB antibody should be determined by the researcher.

Immunogen

Amino acids RQKDPVKGMTADD were used as the immunogen for this PITPNB antibody.

Storage

Aliquot and store the PITPNB antibody at -20oC.