


PITPNM3 Antibody (R36107)

Catalog No.	Formulation	Size
R36107-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	83394
Applications	Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This PITPNM3 antibody is available for research use only.



IHC testing of FFPE human cerebral cortex with PITPNM3 antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: PYK2 N-terminal domain-interacting receptor 1, NIR1

Application Notes

Optimal dilution of the PITPNM3 antibody should be determined by the researcher.

Immunogen

Amino acids SWARGPPKFESVP were used as the immunogen for this PITPNM3 antibody.

Storage

Aliquot and store the PITPNM3 antibody at -20oC.