


PITX2 Antibody (RQ5602)

Catalog No.	Formulation	Size
RQ5602	0.5mg/ml if reconstituted with 0.2ml sterile DI water	100 ug

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit IgG
Purity	Affinity purified
Buffer	Lyophilized from 1X PBS with 2% Trehalose and 0.025% sodium azide
UniProt	Q99697
Applications	Western Blot : 0.25-0.5ug/ml
Limitations	This PITX2 antibody is available for research use only.



Western blot testing of human 1) Caco-2, 2) HEK293, 3) U-2 OS, 4) HeLa, 5) A549, 6) U-87 MG, 7) rat heart and 8) mouse RAW264.7 lysate with PITX2 antibody. Predicted molecular weight ~35 kDa.

Description

Paired-like homeodomain transcription factor 2 also known as pituitary homeobox 2 is a protein that in humans is encoded by the PITX2 gene. It is mapped to 4q25. This gene encodes a member of the RIEG/PITX homeobox family, which is in the bicoid class of homeodomain proteins. The encoded protein acts as a transcription factor and regulates procollagen lysyl hydroxylase gene expression. This protein plays a role in the terminal differentiation of somatotroph and lactotroph cell phenotypes, is involved in the development of the eye, tooth and abdominal organs, and acts as a transcriptional regulator involved in basal and hormone-regulated activity of prolactin. Mutations in this gene are associated with Axenfeld-Rieger syndrome, iridogoniodysgenesis syndrome, and sporadic cases of Peters anomaly. A similar protein in other vertebrates is involved in the determination of left-right asymmetry during development.

Alternatively spliced transcript variants encoding distinct isoforms have been described.

Application Notes

Optimal dilution of the PITX2 antibody should be determined by the researcher.

Immunogen

Amino acids METNCRKLVSACVQLGVQPAAVECLFSKDSEIKK were used as the immunogen for the PITX2 antibody.

Storage

After reconstitution, the PITX2 antibody can be stored for up to one month at 4oC. For long-term, aliquot and store at -20oC. Avoid repeated freezing and thawing.