


PLA2 Antibody / PLA2G2A (R36113)

Catalog No.	Formulation	Size
R36113-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5320
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This PLA2 antibody is available for research use only.



Western blot testing of placenta lysate with PLA2 antibody at 0.1ug/ml. Predicted molecular weight ~16kDa.

Description

Additional name(s) for this target protein: secretory Phospholipase A2; phospholipase A2, group IIA (PLA2G2A)

Application Notes

Optimal dilution of the PLA2 antibody should be determined by the researcher.

Immunogen

Amino acids SYKFSNSGSRIT were used as the immunogen for this PLA2 antibody.

Storage

Aliquot and store the PLA2 antibody at -20oC.