


PLA2R1 Antibody (R33831)

Catalog No.	Formulation	Size
R33831-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	22925
Applications	Flow Cytometry : 10ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This PLA2R1 antibody is available for research use only.



FACS testing of fixed and permeabilized human HEK293 cells with PLA2R1 antibody (blue) at 10ug/ml and naive goat Ig (black).

Description

Additional name(s) for this target protein: Phospholipase A2 receptor 1

Application Notes

Optimal dilution of the PLA2R1 antibody should be determined by the researcher.

Immunogen

Amino acids EKKKDTPKQHGT were used as the immunogen for this PLA2R1 antibody.

Storage

Aliquot and store the PLA2R1 antibody at -20oC.