


PMEL17 Antibody (R33299)

Catalog No.	Formulation	Size
R33299-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	6490
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This PMEL17 antibody is available for research use only.



Western blot testing of human skin lysate with PMEL17 antibody at 1ug/ml. Predicted molecular weight: ~70/100/26 kDa (unmodified/glycosylated/M-beta fragment).

Description

Additional name(s) for this target protein: Melanocyte protein mel 17; SILV, gp100

Application Notes

Optimal dilution of the PMEL17 antibody should be determined by the researcher.

Immunogen

Amino acids CPIGENSPLLSGQQ were used as the immunogen for this PMEL17 antibody.

Storage

Aliquot and store the PMEL17 antibody at -20oC.

References (1)