


Podoplanin Antibody (R33852)

Catalog No.	Formulation	Size
R33852-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	10630
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:1000
Limitations	This Podoplanin antibody is available for research use only.



Western blot testing of human cerebellum lysate with Podoplanin antibody at 1ug/ml.
Predicted molecular weight ~25/36kDa (unmodified/glycosylated).

Description

Additional name(s) for this target protein: PDPN; AGGRUS; Glycoprotein 36 (GP36); GP38; GP40; HT1A-1; OTS8; PA2.26

Application Notes

Optimal dilution of the Podoplanin antibody should be determined by the researcher.

Immunogen

Amino acids EKVDGDTQTTVEKD were used as the immunogen for this Podoplanin antibody.

Storage

Aliquot and store the Podoplanin antibody at -20oC.