


POLG Antibody / DNA polymerase subunit gamma 1 [clone 24P93] (RQ8888)

Catalog No.	Formulation	Size
RQ8888	Antibody in PBS with 0.02% sodium azide, 50% glycerol and 0.4-0.5mg/ml BSA	100 ul

Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	24P93
Purity	Affinity chromatography
UniProt	P54098
Applications	Western Blot : 1:500-1:2000
Limitations	This POLG antibody is available for research use only.



Western blot testing of human 1) MCF7, 2) Jurkat and 3) SiHa cell lysate with POLG antibody. Predicted molecular weight ~140 kDa.