


PON2 Antibody / Paraoxonase 2 (R33967)

Catalog No.	Formulation	Size
R33967-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Rat, Dog, Pig, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5445
Localization	Cytoplasmic
Applications	Immunohistochemistry (FFPE) : 2-4ug/ml ELISA (peptide) LOD : 1:4000
Limitations	This PON2 antibody is available for research use only.



IHC testing of FFPE human liver with PON2 antibody at 2ug/ml. HIER: steamed with pH6 citrate buffer, HRP-staining.

Description

Additional name(s) for this target protein: Paraoxonase 2

Application Notes

Optimal dilution of the PON2 antibody should be determined by the researcher.

1. This PON2 antibody will detect isoforms 1 and 2 of the protein.

Immunogen

Amino acids ERLALRNRLKAS were used as the immunogen for this PON2 antibody.

Storage

Aliquot and store the PON2 antibody at -20oC.