


POR Antibody / CYPOR / Cytochrome P450 reductase (R36231)

Catalog No.	Formulation	Size
R36231-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse, Rat
Predicted Reactivity	Cow, Dog, Pig, Rabbit
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5447
Localization	Cytoplasmic, membrane
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:16000
Limitations	This POR antibody is available for research use only.



Western blot testing of rat liver lysate with POR antibody at 0.1ug/ml. Predicted molecular weight: ~77kDa, observed here at ~90kDa.

Description

Additional name(s) for this target protein: Cytochrome P450 reductase; CYPOR; CPR

Application Notes

Optimal dilution of the POR antibody should be determined by the researcher.

Immunogen

Amino acids TNPPRTNVLYE were used as the immunogen for this POR antibody.

Storage

Aliquot and store the POR antibody at -20oC.