


Pou3f3a Antibody (F52241)

Catalog No.	Formulation	Size
F52241-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F52241-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	P56224
Applications	Western Blot : 1:1000
Limitations	This Pou3f3a antibody is available for research use only.



Description

Transcription factor that may play important roles in patterning the embryonic brain.

Application Notes

Titration of the Pou3f3a antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A portion of amino acids 381-407 from the zebrafish protein was used as the immunogen for this Pou3f3a antibody.

Storage

Aliquot the Pou3f3a antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.