


POU4F3 Antibody (R33696)

Catalog No.	Formulation	Size
R33696-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Dog, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5459
Applications	Western Blot : 0.2-0.6ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This POU4F3 antibody is available for research use only.



Western blot testing of human heart lysate with POU4F3 antibody at 0.2ug/ml. Predicted molecular weight: ~37kDa.

Description

Additional name(s) for this target protein: POU class 4 homeobox 3

Application Notes

Optimal dilution of the POU4F3 antibody should be determined by the researcher.

1. This POU4F3 antibody does not cross react with similar proteins POU4F1 and POU4F2.

Immunogen

Amino acids EAAYREKNSKPE were used as the immunogen for this POU4F3 antibody.

Storage

Aliquot and store the POU4F3 antibody at -20oC.