


## PRAM1 Antibody (F53937)

Catalog No.	Formulation	Size
F53937-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

**Bulk quote request**

<b>Availability</b>	1-3 business days
<b>Species Reactivity</b>	Human
<b>Format</b>	Antigen affinity purified
<b>Clonality</b>	Polyclonal (rabbit origin)
<b>Isotype</b>	Rabbit Ig
<b>Purity</b>	Antigen affinity
<b>UniProt</b>	Q96QH2
<b>Applications</b>	Western Blot : 1:500-1000
<b>Limitations</b>	This PRAM1 antibody is available for research use only.



Western blot testing of human THP-1 cell lysate with PRAM1 antibody at 1:1000.  
Predicted molecular weight: 79 kDa but can also be observed at ~97 kDa.

### Description

PML-RARA-regulated adapter molecule 1 may be involved in myeloid differentiation. May be involved in integrin signaling in neutrophils. Binds to PtdIns(4)P.

### Application Notes

Titration of the PRAM1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

### Immunogen

A portion of amino acids 587-618 from the human protein was used as the immunogen for the PRAM1 antibody.

### Storage

Aliquot the PRAM1 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.