


PRAM1 Antibody (R33436)

Catalog No.	Formulation	Size
R33436-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Pig
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	84106
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This PRAM1 antibody is available for research use only.



Western blot testing of Jurkat lysate with PRAM1 antibody at 2ug/ml. Predicted molecular weight: ~79kDa but routinely observed at 97kDa.

Description

Additional name(s) for this target protein: PML-RARA target gene encoding an Adaptor Molecule-1

Application Notes

Optimal dilution of the PRAM1 antibody should be determined by the researcher.

Immunogen

Amino acids DFCDPLENQPLPLGR were used as the immunogen for this PRAM1 antibody.

Storage

Aliquot and store the PRAM1 antibody at -20oC.