


PRDM16 Antibody [clone 1740CT793.52.44] (F53945)

Catalog No.	Formulation	Size
F53945-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F53945-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG1, kappa
Clone Name	1740CT793.52.44
Purity	Protein G affinity chromatography
UniProt	Q9HAZ2
Applications	Western Blot : 1:2000-4000
Limitations	This PRDM16 antibody is available for research use only.



Western blot testing of human 1) Jurkat and 2) K562 cell lysate with PRDM16 antibody at 1:4000. Predicted molecular weight: 140 kDa.

Description

PR domain zinc finger protein 16 binds DNA and functions as a transcriptional regulator. Functions in the differentiation of brown adipose tissue (BAT) which is specialized in dissipating chemical energy in the form of heat in response to cold or excess feeding while white adipose tissue (WAT) is specialized in the storage of excess energy and the control of systemic metabolism. Together with CEBPB, regulates the differentiation of myoblastic precursors into brown adipose cells. Functions also as a repressor of TGF-beta signaling. Isoform 4 may regulate granulocytes differentiation.

Application Notes

Titration of the PRDM16 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A partial recombinant protein corresponding to amino acids 779-996 from the human protein was used as the immunogen for the PRDM16 antibody.

Storage

Aliquot the PRDM16 antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.