


PRKCDBP Antibody (R35755)


Catalog No.	Formulation	Size
R35755-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Cow, Dog, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	112464
Applications	Western Blot : 0.1-0.3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This PRKCDBP antibody is available for research use only.



Western blot testing of human adipose cell lysate with PRKCDBP antibody at 0.1ug/ml. Predicted molecular weight: ~28kDa, routinely observed at 35-40kDa.



Western blot testing of mouse ovary (A) and rat uterus (B) lysate with PRKCDBP antibody at 0.3ug/ml. Predicted molecular weight: ~28kDa, routinely observed at 35-40kDa.

Description

Additional name(s) for this target protein: Protein kinase C delta binding protein

Application Notes

Optimal dilution of the PRKCDBP antibody should be determined by the researcher.

Immunogen

Amino acids RRTGLQKVQSLR were used as the immunogen for this PRKCDBP antibody.

Storage

Aliquot and store the PRKCDBP antibody at -20oC.