

PRMT10 Antibody / PRMT9 / Protein arginine N-methyltransferase 9 (F46110)

Catalog No.	Formulation	Size
F46110-0.4ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.4 ml
F46110-0.08ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.08 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human, Mouse
Format	Antigen affinity purified
Clonality	Polyclonal (rabbit origin)
Isotype	Rabbit Ig
Purity	Antigen affinity
UniProt	Q6P2P2
Applications	Western Blot : 1:1000
Limitations	This PRMT10 antibody is available for research use only.

250
130
95
72
55

Western blot testing of mouse cerebellum tissue lysate with PRMT10 antibody. Predicted molecular weight ~95 kDa.

250
130
95
72
55

Western blot testing of human Jurkat cell lysate with PRMT10 antibody. Predicted molecular weight ~95 kDa.

