


## PRMT7 Antibody (R35135)

Catalog No.	Formulation	Size
R35135-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Mouse
Predicted Reactivity	Dog, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	54496
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:128000
Limitations	This PRMT7 antibody is available for research use only.



Western blot testing of HeLa lysate with PRMT7 antibody at 0.3ug/ml. Predicted molecular weight: ~78kDa.

## Description

Additional name(s) for this target protein: Protein arginine methyltransferase 7

## Application Notes

Optimal dilution of the PRMT7 antibody should be determined by the researcher.

## Immunogen

Amino acids PRFGEINDQDRTDR were used as the immunogen for this PRMT7 antibody.

## Storage

Aliquot and store the PRMT7 antibody at -20oC.