


Prokineticin receptor 2 Antibody (R34807)


Catalog No.	Formulation	Size
R34807-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human, Rat
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	128674
Applications	Western Blot : 0.3-1ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This Prokineticin receptor 2 antibody is available for research use only.



Western blot of human cerebellum (A), rat brain (B) and rat testis (C) lysate with Prokineticin receptor 2 antibody at 0.1ug/ml. Predicted molecular weight (human): ~44kDa.



Western blot testing of rat stomach lysate with Prokineticin receptor 2 antibody at 1ug/ml. Predicted molecular weight (rat): ~47kDa.

Description

Additional name(s) for this target protein: PROKR2

Application Notes

Optimal dilution of the Prokineticin receptor 2 antibody should be determined by the researcher.

Immunogen

Amino acids NGNTSFTPNFN were used as the immunogen for this Prokineticin receptor 2 antibody.

Storage

Aliquot and store the Prokineticin receptor 2 antibody at -20oC.