


PSMB10 Antibody (R35169)


Catalog No.	Formulation	Size
R35169-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5699
Applications	Western Blot : 0.3-1ug/ml Immunohistochemistry (FFPE) : 5-10ug/ml ELISA (peptide) LOD : 1:8000
Limitations	This PSMB10 antibody is available for research use only.



Western blot testing of human lung lysate with PSMB10 antibody at 0.3ug/ml. Predicted molecular weight: ~29kDa.



IHC staining of FFPE human tonsil with PSMB10 antibody at 5ug/ml. HIER: steamed with pH6 citrate buffer, AP-staining.

Description

Additional name(s) for this target protein: Proteasome (prosome, macropain) subunit, beta type, 10

Application Notes

Optimal dilution of the PSMB10 antibody should be determined by the researcher.

Immunogen

Amino acids PTEPVKRSGRYH were used as the immunogen for this PSMB10 antibody.

Storage

Aliquot and store the PSMB10 antibody at -20oC.