


PTF1A Antibody (R36349)

Catalog No.	Formulation	Size
R36349-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Dog, Mouse, Rat, Zebrafish
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	256297
Applications	Western Blot : 0.5-1.5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This PTF1A antibody is available for research use only.



Western blot testing of human cerebellum lysate with PTF1A antibody at 0.5ug/ml.
Predicted molecular weight: ~35kDa, routinely observed at ~48kDa. (Ref 1)

Description

Additional name(s) for this target protein: Pancreas specific transcription factor 1a; PTF1-p48; bHLH transcription factor p48

Application Notes

Optimal dilution of the PTF1A antibody should be determined by the researcher.

Immunogen

Amino acids WTDEKQLKEQN were used as the immunogen for this PTF1A antibody.

Storage

Aliquot and store the PTF1A antibody at -20oC.

References (1)