


## PUM2 Antibody (R33599)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R33599-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                      |                                                           |
|----------------------|-----------------------------------------------------------|
| Availability         | 1-3 business days                                         |
| Species Reactivity   | Human                                                     |
| Predicted Reactivity | Cow                                                       |
| Format               | Antigen affinity purified                                 |
| Clonality            | Polyclonal (goat origin)                                  |
| Isotype              | Goat Ig                                                   |
| Purity               | Antigen affinity                                          |
| Gene ID              | 23369                                                     |
| Applications         | Western Blot : 0.5-2ug/ml<br>ELISA (peptide) LOD : 1:4000 |
| Limitations          | This PUM2 antibody is available for research use only.    |



Western blot testing of HEK293 lysate overexpressing human PUM2-FLAG with PUM2 antibody (0.5ug/ml) in Lane A and probed with anti-FLAG (1/3000) in lane C. Mock-transfected HEK293 probed with PUM2 antibody (1ug/ml) in Lane B.

## Description

Additional name(s) for this target protein: Pumilio homolog 2

## Application Notes

Optimal dilution of the PUM2 antibody should be determined by the researcher.

## Immunogen

Amino acids DQNRDLKQGDDDDSK were used as the immunogen for this PUM2 antibody.

## Storage

Aliquot and store the PUM2 antibody at -20oC.