


## PYROXD1 Antibody (R36135)

| Catalog No.  | Formulation                                                                   | Size   |
|--------------|-------------------------------------------------------------------------------|--------|
| R36135-100UG | 0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide | 100 ug |

[Bulk quote request](#)

|                    |                                                           |
|--------------------|-----------------------------------------------------------|
| Availability       | 1-3 business days                                         |
| Species Reactivity | Human                                                     |
| Format             | Antigen affinity purified                                 |
| Clonality          | Polyclonal (goat origin)                                  |
| Isotype            | Goat Ig                                                   |
| Purity             | Antigen affinity                                          |
| Gene ID            | 79912                                                     |
| Applications       | Western Blot : 1-3ug/ml<br>ELISA (peptide) LOD : 1:128000 |
| Limitations        | This PYROXD1 antibody is available for research use only. |



Western blot testing of human colorectal cancer lysate with PYROXD1 antibody at 1ug/ml. Predicted/observed molecular weight ~56kDa. The ~48kDa band is likely a breakdown product.

### Description

Additional name(s) for this target protein: Pyridine nucleotide-disulphide oxidoreductase domain 1; Pyridine nucleotide-disulfide oxidoreductase domain-containing protein 1

### Application Notes

Optimal dilution of the PYROXD1 antibody should be determined by the researcher.

### Immunogen

Amino acids TCAEQLATHFPSED were used as the immunogen for this PYROXD1 antibody.

## Storage

Aliquot and store the PYROXD1 antibody at -20oC.