


RAB5C Antibody [clone 1616CT314.65.31] (F53802)


Catalog No.	Formulation	Size
F53802-0.2ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.2 ml
F53802-0.05ML	In 1X PBS, pH 7.4, with 0.09% sodium azide	0.05 ml

Bulk quote request

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Monoclonal (mouse origin)
Isotype	Mouse IgG2b, kappa
Clone Name	1616CT314.65.31
Purity	Protein G affinity chromatography
UniProt	P51148
Applications	Western Blot : 1:1000-2000 FACS (intracellular) : 1:25
Limitations	This RAB5C antibody is available for research use only.



Western blot testing of human 1) A431, 2) HL-60 and 3) U2OS cell lysate with RAB5C antibody at 1:2000. Predicted molecular weight: 23 kDa.



Intracellular staining of U2OS cells with RAB5C antibody (green) and [isotype control](#) (blue) at 1:25 dilution.

Description

RAB5C is involved in protein transport and likely involved in vesicular traffic.

Application Notes

Titration of the RAB5C antibody may be required due to differences in protocols and substrate sensitivity.

Immunogen

Human recombinant protein was used as the immunogen for this RAB5C antibody.

Storage

Aliquot the RAB5C antibody and store frozen at -20oC or colder. Avoid repeated freeze-thaw cycles.