

Rabenosyn 5 Antibody (R33864)

Catalog No.	Formulation	Size
R33864-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	64145
Applications	Western Blot : 0.3-1.5 ug/ml ELISA (peptide) LOD : 1:64000
Limitations	This Rabenosyn 5 antibody is available for research use only.



Western blot testing of A431 lysate with Rabenosyn 5 antibody at 0.3ug/ml. Predicted molecular weight: ~89kDa but commonly observed at ~110kDa.

Description

Additional name(s) for this target protein: FYVE-finger-containing Rab5 effector protein rabenosyn-5; zinc finger FYVE domain containing 20; ZFYVE20

Application Notes

Optimal dilution of the Rabenosyn 5 antibody should be determined by the researcher.

Immunogen

Amino acids ELKHTLAKQKGGTD were used as the immunogen for this Rabenosyn 5 antibody.

Storage

Aliquot and store the Rabenosyn 5 antibody at -20oC.