


## RAD51C Antibody (R35618)

Catalog No.	Formulation	Size
R35618-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	5889
Applications	Western Blot : 1-3ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This RAD51C antibody is available for research use only.



Western blot testing of HeLa lysate with RAD51C antibody at 2ug/ml. Predicted molecular weight: ~42kDa.

## Description

Additional name(s) for this target protein: RAD51 homolog C, RAD51 homolog 3

## Application Notes

Optimal dilution of the RAD51C antibody should be determined by the researcher.

## Immunogen

Amino acids RGKTFRFEMQRDL were used as the immunogen for this RAD51C antibody.

## Storage

Aliquot and store the RAD51C antibody at -20oC.