


RANBP16 Antibody (R34927)


Catalog No.	Formulation	Size
R34927-100UG	0.5 mg/ml in 1X TBS, pH7.3, with 0.5% BSA (US sourced) and 0.02% sodium azide	100 ug

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Predicted Reactivity	Mouse, Dog, Pig, Cow
Format	Antigen affinity purified
Clonality	Polyclonal (goat origin)
Isotype	Goat Ig
Purity	Antigen affinity
Gene ID	23039
Localization	Cytoplasmic, nuclear
Applications	Western Blot : 0.5-1ug/ml Immunohistochemistry (FFPE) : 5ug/ml ELISA (peptide) LOD : 1:32000
Limitations	This RANBP16 antibody is available for research use only.



Western blot testing of HeLa lysate with RANBP16 antibody at 0.5ug/ml. Predicted molecular weight: ~124 kDa.



IHC testing of FFPE human testis tissue with RANBP16 antibody at 5ug/ml. Required HIER: steamed antigen retrieval with pH6 citrate buffer; AP-staining.

Description

Additional name(s) for this target protein: Exportin 7; XPO7

Application Notes

Optimal dilution of the RANBP16 antibody should be determined by the researcher.

Immunogen

Amino acids NSTYGVNSNDMMS were used as the immunogen for this RANBP16 antibody.

Storage

Aliquot and store the RANBP16 antibody at -20oC.