


Recombinant Adipophilin Antibody / ADRP / Perilipin 2 / PLIN2 [clone ADFP/2755R] (V3970)

Catalog No.	Formulation	Size
V3970-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	100 ug
V3970-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide	20 ug
V3970SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug
V3970IHC-7ML	Prediluted in 1X PBS with 0.1 mg/ml BSA (US sourced) and 0.05% sodium azide; *For IHC use only*	7 ml


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)


Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG, kappa
Clone Name	ADFP/2755R
Purity	Protein A affinity chromatography
UniProt	Q99541
Localization	Cytoplasmic
Applications	ELISA (order BSA/sodium Azide-free Format For Coating) : Western Blot : 1-2ug/ml Immunohistochemistry (FFPE) : 1-2ug/ml for 30 min at RT
Limitations	This recombinant Adipophilin antibody is available for research use only.



IHC testing of FFPE human adrenal tissue with recombinant Adipophilin antibody (clone ADFP/2755R). Required HIER: steam sections in pH6, 10mM citrate buffer for 10-20 min.



Western blot testing of human liver lysate with recombinant Adipophilin antibody.
Predicted molecular weight ~48 kDa.



SDS-PAGE analysis of purified, BSA-free recombinant Adipophilin antibody (clone ADFP/2755R) as confirmation of integrity and purity.

Description

Recognizes a protein of 48kDa, which is identified as Adipophilin. It belongs to the perilipin family, members of which coat intracellular lipid storage droplets. This protein is associated with the lipid globule surface membrane material, and maybe involved in development and maintenance of adipose tissue. However, it is not restricted to adipocytes as previously thought, but is found in a wide range of cultured cell lines, including fibroblasts, endothelial and epithelial cells, and tissues, such as lactating mammary gland, adrenal cortex, Sertoli and Leydig cells, and hepatocytes in alcoholic liver cirrhosis, suggesting that it may serve as a marker of lipid accumulation in diverse cell types and diseases.

Application Notes

Titering of the recombinant Adipophilin antibody may be required for optimal performance.

1. The prediluted format is supplied in a dropper bottle and is optimized for use in IHC. After epitope retrieval step (if required), drip mAb solution onto the tissue section and incubate at RT for 30 min.

Immunogen

A portion of amin acids 249-376 from the human protein were used as the immunogen for this recombinant Adipophilin antibody.

Storage

The recombinant Adipophilin antibody (with azide) can be stored at 2-8oC. The azide-free format should be aliquoted and stored at -20oC or colder.