


Recombinant Arginase 1 Antibody [clone RM377] (R20396)

Catalog No.	Formulation	Size
R20396-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM377
Purity	Protein A purified from animal origin-free supernatant
UniProt	P05089
Localization	Cytoplasmic
Applications	Immunohistochemistry (FFPE) : 1:1000-1:2000 Western Blot : 1:5000-1:10000
Limitations	This recombinant Arginase 1 antibody is available for research use only.



IHC staining of FFPE human liver cancer tissue with recombinant Arginase-1 antibody at 1:1250.



Western blot testing of human liver tissue with recombinant Arginase 1 antibody at 1:10,000 dilution. Predicted molecular weight ~35 kDa.

Description

This antibody reacts to human Arginase-1 (ARG1).

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant Arginase 1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A peptide corresponding to the C-terminus of human Arginase-1 (ARG1) was used as the immunogen for the recombinant Arginase 1 antibody.

Storage

Store the recombinant Arginase-1 antibody at -20oC.