

Recombinant BAX Antibody [clone BAX/3158R] (V5783)

Catalog No.	Formulation	Size
V5783-100UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	100 ug
V5783-20UG	0.2 mg/ml in 1X PBS with 0.1 mg/ml BSA (US sourced), 0.05% sodium azide	20 ug
V5783SAF-100UG	1 mg/ml in 1X PBS; BSA free, sodium azide free	100 ug


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)


Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG, kappa
Clone Name	BAX/3158R
UniProt	Q07812
Localization	Cytoplasm
Applications	Western Blot : 2-4ug/ml
Limitations	This recombinant BAX antibody is available for research use only.



Western blot testing of human A549 and MCF7 cell lysate with recombinant BAX antibody (clone BAX/3158R). Predicted molecular weight ~21 kDa.



Western blot testing of human A549 cell lysate with recombinant BAX antibody (clone BAX/3158R). Predicted molecular weight ~21 kDa.



SDS-PAGE analysis of purified, BSA-free recombinant BAX antibody (clone BAX/3158R) as confirmation of integrity and purity.

Description

Recognizes a protein of 21 kDa, identified as the Bax protein. This MAb is highly specific to Bax and shows no cross-reaction with Bcl-2 or Bcl-X protein. Bcl-2 blocks cell death following a variety of stimuli. Bax has extensive amino acid homology with Bcl-2 and it homodimerizes and forms heterodimers with Bcl-2. Overexpression of Bax accelerates apoptotic death induced by cytokine deprivation in an IL-3 dependent cell line, and Bax also counters the death repressor activity of Bcl-2.

Application Notes

Optimal dilution of the recombinant BAX antibody should be determined by the researcher.

Immunogen

A portion of amino acids 1-100 from human BAX protein was used as the immunogen for the recombinant BAX antibody.

Storage

Aliquot the recombinant BAX antibody and store frozen at -20°C or colder. Avoid repeated freeze-thaw cycles.