


Recombinant BOB-1 Antibody / POU2AF1 [clone RM378] (R20397)

Catalog No.	Formulation	Size
R20397-0.1ML	Antibody in PBS with 50% glycerol, 1% BSA and 0.09% sodium azide	100 ul


Recombinant **RABBIT MONOCLONAL**

[Bulk quote request](#)

Availability	1-3 business days
Species Reactivity	Human
Format	Purified
Clonality	Recombinant Rabbit Monoclonal
Isotype	Rabbit IgG
Clone Name	RM378
Purity	Protein A purified from animal origin-free supernatant
UniProt	Q16633
Localization	Nuclear
Applications	Immunohistochemistry (FFPE) : 1:1000-1:2000 Western Blot : 1:5000-1:10000
Limitations	This recombinant BOB-1 antibody is available for research use only.



Western blot testing of human Daudi cell lysate with recombinant BOB-1 antibody at 1:10,000 dilution. Predicted molecular weight: ~28 kDa (unmodified), 35-40 kDa (ubiquitinated).



IHC staining of FFPE human tonsil tissue with recombinant BOB-1 antibody at 1:1250.

Description

This antibody reacts to human BOB-1/OBF-1 (B-cell-specific coactivator) (POU domain class 2-associating factor 1).

Application Notes

The stated application concentrations are suggested starting points. Titration of the recombinant BOB-1 antibody may be required due to differences in protocols and secondary/substrate sensitivity.

Immunogen

A peptide corresponding to the C-terminus of human BOB-1/OBF-1 was used as the immunogen for the recombinant BOB-1 antibody.

Storage

Store the recombinant BOB-1 antibody at -20oC.